

Лабораторные расходные материалы и посуда

Описание

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: azp@nt-rt.ru || сайт: <https://labor-pilz.nt-rt.ru/>

Additional core for adhesive tape dispenser

Order-No.: 01175



approx. 25 mm (1 inch) / approx. 76 mm (3 inches)

With additional cores, more adhesive tapes or labeling bands can be installed in the donor for adhesive tapes at the same time. Each core enables the installation of an adhesive tape or lettering band with a core hole diameter of approx. 25 mm (1 inch) or approx. 76 mm (3 inches).

Compliance information

WEEE obligation..... No

Slotted rail mini with handles

Order-No.: 01602



The Nutschiene can be used in many positions in Smartrack®. A larger number of Nutschienen offers more options for the placement of modules through an optimal vertical arrangement and thus offers more scope for personalizing the rack and adapting the module arrangement to the individual needs of each individual user.

Compliance information

WEEE obligation..... No

Slotted rail with handles

Order-No.: 01601



The Nutschiene can be used in many positions in Smartrack®. A larger number of Nutschienen offers more options for the placement of modules through an optimal vertical arrangement and thus offers more scope for personalizing the rack and adapting the module arrangement to the individual needs of each individual user.

Compliance information

WEEE obligation..... No

Chemistry module set

Order-No.: 01503



The chemistry module set was specially designed for work in the field of chemistry and includes the following components (article number of the individual module for better understanding in brackets):

1 x Nutschiene with handles (Art.No. 01601)

9 modules:

- 3 x bracket for chromatography pillars (Art.No. 01107)
- 3 x bracket for glass appliances with NS29 (Art.No. 01114)
- 1 x bracket for glass appliances with NS29, extended (Art.No. 01116)
- 2 x bracket for glass appliances with NS14 (Art.No.01111)

Compliance information

WEEE obligation..... No

Biology module set

Order-No.: 01502



The module set biology was specially designed for work in the field of biology and includes the following components (article number of the individual module for better understanding in brackets):

1 x Nutschiene with handles (Art.No. 01601)

8 modules:

- 1 x shelf (Art.No. 01142)

- 1 x shelter insert- for centrifuge tube (Art.No. 01102)

- 2 x bracket for laboratory thread bottles 100 - 1000 ml (GL45) (Art.No. 01122)

- 1 x bracket for 3 pipettes with hooks (Art.No. 01145)

- 3 x bracket for pipette tip boxes (Art.No. 01147)

Compliance information

WEEE obligation..... No

Basic equipment module set

Order-No.: 01501



The basic equipment module set is a basic equipment for the organization of basic accessories in everyday laboratory life and includes the following components (article number of the individual module for better understanding in brackets):

6 modules:

- 1 x shelf (Art.No. 01142)
- 1 x shelf use- for small accessories (item no. 01148)
- 1 x bracket for garbage bags (Art.No. 01144)
- 1 x bracket for glove / tissue boxes (Art.No. 01143)
- 1 x bracket for spray bottles (Art.No. 01146)
- 1 x wall bracket (Art.No. 01126)

Compliance information

WEEE obligation..... No

Magnet for pipette tip boxes

Order-No.: 01151



The magnet can be inserted into a pipette box and thus enables the pipette tip box on the bracket for pipette tip boxes.

Compliance information

WEEE obligation..... No

Holder for pipetting aid

Order-No.: 01178



Type: For Brand Accu-Jet® Pro

The module enables the safe hanging of a pipetting help.

Compliance information

WEEE obligation..... No

Holder for chromatography columns

Order-No.: 01177



Type: Ø 5 - 10 mm

The module enables an effortless hanging of chromatography pillars. The bracket adapts to the scope of columns through a snap mechanism and thus offers a safe hold. From a length of more than 15 cm of the hung device, it is advisable to use two modules arranged one above the other to ensure the stability of the structure.

Compliance information

WEEE obligation..... No

Holder for chromatography columns

Order-No.: 01176



Type: Ø 10 - 15 mm

The module enables an effortless hanging of chromatography pillars. The bracket adapts to the scope of columns through a snap mechanism and thus offers a safe hold. From a length of more than 15 cm of the hung device, it is advisable to use two modules arranged one above the other to ensure the stability of the structure.

Compliance information

WEEE obligation..... No

Holder for needle ejector container

Order-No.: 01174



For Vernacare Sharpsafe® 1L

The module enables easy hits of a cannula drop box.

Compliance information

WEEE obligation..... No

Holder for syringe bottles

Order-No.: 01173



For 1000 ml (up to Ø 95 mm)

The module enables easy hanging of spray bottles. The module can also be used to store screw cap bottles of similar sizes.

Compliance information

WEEE obligation..... No

Holder for pipettes

Order-No.: 01168



For Eppendorf Reference® 2

The module enables an easy hanging of an Eppendorf Reference® 2 pipette. The pipette tips are not loaded and therefore not contaminated.

Compliance information

WEEE obligation..... No

Holder for pipettes

Order-No.: 01166



For Gilson® Pipetman® Classic / G

The module enables a simple hanging of a Gilson® Pipetman® Classic or Gilson® Pipetman® G Pipette. The pipette tips are not loaded and therefore not contaminated.

Compliance information

WEEE obligation..... No

Holder for chromatography columns

Order-No.: 01165



Type: Ø 15 - 25 mm

The module enables an effortless hanging of chromatography pillars. The bracket adapts to the scope of columns through a snap mechanism and thus offers a safe hold. From a length of more than 15 cm of the hung device, it is advisable to use two modules arranged one above the other to ensure the stability of the structure.

Compliance information

WEEE obligation..... No

Holder for pipettes

Order-No.: 01164



For Eppendorf Reference®

The module enables an easy hanging of an Eppendorf® Reference Pipette. The pipette tips are not loaded and therefore not contaminated.

Compliance information

WEEE obligation..... No

Holder for barcode scanners

Order-No.: 01163



Type: i.a. Suitable for Newland NLS-HR22 and Aplic 304674

The module enables a simply hanging of a barcode scanner.

Compliance information

WEEE obligation..... No

Holder for kitchen rolls

Order-No.: 01162



The module enables easy hanging of kitchen rolls. The roll can rotate in the module, which enables easy rolling of individual or more leaves.

Compliance information

WEEE obligation..... No

Holder for open fronted storage bins

Order-No.: 01158



Type: i.a. Suitable for Schäfer LF110, LF210, LF211, LF220

The module enables the safe hanging of visual storage boxes or stacking boxes with a hanging edge, which typically serves as a holding device for wall rails and shelf systems.

Compliance information

WEEE obligation..... No

Holder for pipetting aid

Order-No.: 01157



Type: Brand Accu-Jet® S

The module enables the safe hanging of a pipetting help.

Compliance information

WEEE obligation..... No

Holder for pipetting aid

Order-No.: 01156



Type: Integra-Pipetboy Pro

The module enables the safe hanging of a pipetting help.

Compliance information

WEEE obligation..... No

Holder for microliter syringes

Order-No.: 01155



The module enables the safe hanging of three microliter injections. Due to the snap mechanism, the hung microliter injections have a safe hold. The shape of the module shields the cannulas of the microliter injections and the user thus protected against unintentional touch with the top of the cannula.

Compliance information

WEEE obligation..... No

Holder for pipetting aid

Order-No.: 01154



Type: Integra pipetboy 2

The module enables the safe hanging of a pipetting help.

Compliance information

WEEE obligation..... No

Holder for separatory funnel / shaking funnel

Order-No.: 01153



Type: Ø> 105 mm

The module enables a simple hanging of vendors / shüttstelter.

A notice:

The diameter of the divorce is measured at its broadest place

Compliance information

WEEE obligation..... No

Holder for separating funnel / shaking funnel

Order-No.: 01152



Type: Ø 45 - 85 mm

The module enables a simple hanging of vendors / shüttstelter.

A notice:

The diameter of the divorce is measured at its broadest place

Compliance information

WEEE obligation..... No

Holder for separatory funnel / shaking funnel

Order-No.: 01149



Type: Ø 85 - 105 mm

The module enables a simple hanging of vendors / shüttstelter.

A notice:

The diameter of the divorce is measured at its broadest place

Compliance information

WEEE obligation..... No

Holder for pipette tip boxes

Order-No.: 01147



The module enables easy fixing of pipette tip boxes of different manufacturers and various sizes. A safe hold of the box is guaranteed by a magnet integrated in the module and a magnetic stainless steelron, which is inserted into the pipette tip box. Further stainless steelron can be ordered to work with several pipette tip boxes.

Compliance information

WEEE obligation..... No

Holder for spray bottles

Order-No.: 01146



for 500 ml (up to Ø 75 mm)

The module enables easy hanging of spray bottles. The module can also be used to store screw cap bottles of similar sizes.

Compliance information

WEEE obligation..... No

Holder for pipettes

Order-No.: 01145



For 3 pipets with hooks

The module enables easy hanging of up to three pipets with hooks. The pipette tips are not loaded and therefore not contaminated.

Compliance information

WEEE obligation..... No

Holder for waste bags

Order-No.: 01144



The module enables easy hanging of garbage bags for small waste. Continuation can be prevented by making waste in the immediate area of work.

Compliance information

WEEE obligation..... No

Holder for glove/tissue boxes

Order-No.: 01143



The module enables easy hanging of glove or tissue boxes. The width of the bracket can be set individually and adapted to the size of different boxes. By moving out of the consumables down, this can be easily removed even when the box is low.

Compliance information

WEEE obligation..... No

Holder for Erlenmeyer flasks

Order-No.: 01140



Type: Wide -necked meat cobbles 250 - 1000 ml or corner of the necklace, 2000 - 5000 ml

The module enables easy hanging of Erlenmeyerkroben with a stock edge. Due to the curved front, the Erlenmeyer flock is safely fixed in the module and is protected against falling out.

Compliance information

WEEE obligation..... No

Holder for chromatography columns

Order-No.: 01139



Type: Ø 35 - 46 mm

The module enables an effortless hanging of chromatography pillars. The bracket adapts to the scope of columns through a snap mechanism and thus offers a safe hold. From a length of more than 15 cm of the hung device, it is advisable to use two modules arranged one above the other to ensure the stability of the structure.

Compliance information

WEEE obligation..... No

Holder for well plates

Order-No.: 01125



Well Plates can be inserted into the module and also stacked. The horizontal orientation ensures optimal storage, for example by cell cultures.

Compliance information

WEEE obligation..... No

Holder for tablet

Order-No.: 01124



The module enables simply hanging tablets in portrait and landscape format. The inclination ensures a good view and good usability of the display. The module can be adapted to different tablet sizes via a snap-in mechanism, which can be gradually set.

Compliance information

WEEE obligation..... No

Holder for laboratory flasks

Order-No.: 01122



Type: 100 - 1000ml (GL45)

The module enables easy hanging of laboratory thread bottles. The module is equipped with a rotating cladding ring through which the laboratory thread bottles are also protected against falling out.

Compliance information

WEEE obligation..... No

Holder for laboratory flasks

Order-No.: 01121



Type: 50 ml (GL32)

The module enables easy hanging of laboratory thread bottles. The module is equipped with a rotating cladding ring through which the laboratory thread bottles are also protected against falling out.

Compliance information

WEEE obligation..... No

Holder for laboratory threaded bottles

Order-No.: 01120



Type: 10 - 25 ml (GL25)

The module enables easy hanging of laboratory thread bottles. The module is equipped with a rotating cladding ring through which the laboratory thread bottles are also protected against falling out.

Compliance information

WEEE obligation..... No

Holder for glass appliances with NS29, extended

Order-No.: 01116



The module corresponds to the bracket for glass devices with NS29 in the construction and function. However, the extended variant is particularly suitable for the simplified hanging of e.g. particularly large circular pistons or divorce.

Note on lock tolerances for NS29:

Maximum outer diameter at the narrowest point of the cone cut: 31.4 mm (larger diameter are too wide for the bracket and cannot be absorbed)

Minimal outside diameter at the widest point of the cone cut: 32.8 mm (smaller diameter is too narrow for the bracket and would fall out)

Compliance information

WEEE obligation..... No

Holder for glassware with NS29, plus lock

Order-No.: 01115



The module corresponds to the bracket for glass devices with NS29, but is also equipped with a rotating cladding ring through which the glass devices with standard grinding NS29 are also protected against falling out.

Note on lock tolerances for NS29:

Maximum outer diameter at the narrowest point of the cone cut: 31.4 mm (larger diameter are too wide for the bracket and cannot be absorbed)

Minimal outside diameter at the widest point of the cone cut: 32.8 mm (smaller diameter is too narrow for the bracket and would fall out)

Compliance information

WEEE obligation..... No

Holder for glassware with NS29

Order-No.: 01114



The module enables easy hanging of glass devices with standard grinding NS29. Due to the specially shaped recess, the glass devices are secured against falling out.

Note on lock tolerances for NS29:

Maximum outer diameter at the narrowest point of the cone cut: 31.4 mm (larger diameter are too wide for the bracket and cannot be absorbed)

Minimal outside diameter at the widest point of the cone cut: 32.8 mm (smaller diameter is too narrow for the bracket and would fall out)

Compliance information

WEEE obligation..... No

Holder for glassware with NS14, extended

Order-No.: 01113



The module corresponds to the bracket for glass devices with NS14. However, the extended variant is particularly suitable for the simplified hanging of e.g. particularly large circular pistons or divorce. Note on locklerances for NS14:

Maximum outer diameter at the narrowest point of the cone cut: 17.2 mm (larger diameter are too wide for the bracket and cannot be absorbed)

Minimal outside diameter at the widest point of the cone cut: 18.0 mm (smaller diameter is too narrow for the bracket and would fall out)

Compliance information

WEEE obligation..... No

Holder for glassware with NS14, plus lock

Order-No.: 01112



The module corresponds to the bracket for glass devices with NS14, but is also equipped with a rotating cladding ring through which the glass devices with standard grinding NS14 are also protected against falling out.

Note on locklerances for NS14:

Maximum outer diameter at the narrowest point of the cone cut: 17.2 mm (larger diameter are too wide for the bracket and cannot be absorbed)

Minimal outside diameter at the widest point of the cone cut: 18.0 mm (smaller diameter is too narrow for the bracket and would fall out)

Compliance information

WEEE obligation..... No

Holder for glassware with NS14

Order-No.: 01111



The module enables easy hanging of glass devices with standard grinding NS14. Due to the specially shaped recess, the glass devices are secured against falling out.

Note on locklerances for NS14:

Maximum outer diameter at the narrowest point of the cone cut: 17.2 mm (larger diameter are too wide for the bracket and cannot be absorbed)

Minimal outside diameter at the widest point of the cone cut: 18.0 mm (smaller diameter is too narrow for the bracket and would fall out)

Compliance information

WEEE obligation..... No

Holder for Erlenmeyer flasks

Order-No.: 01108



Type: Wide -necked meat cobbles 25 - 100 ml or corner of the neck of 250 - 500 ml

The module enables easy hanging of Erlenmeyenkroben with a stock edge. Due to the curved front, the Erlenmeyer flock is safely fixed in the module and is protected against falling out.

Compliance information

WEEE obligation..... No

Holder for chromatography columns

Order-No.: 01107



Type: Ø 25 - 35 mm

The module enables an effortless hanging of chromatography pillars. The bracket adapts to the scope of columns through a snap mechanism and thus offers a safe hold. From a length of more than 15 cm of the hung device, it is advisable to use two modules arranged one above the other to ensure the stability of the structure.

Compliance information

WEEE obligation..... No

Fire department experiment case

Order-No.: 104744800



Laboratory mushroom fire brigade experiment-the basis for the combustion and effect of various extinguishing agents of these experimental cases was developed by the fire brigade to convey the basics of combustion and the effect of different extinguishing agents. It is specially intended for use in the case of occupational safety instructions in companies in the laboratory, technology and production sector as support with practical experiments. In addition, the suitcase is used for training by many youth fire brigades. Various experiments were selected for the set, all of which can be carried out in class or training rooms. The attempts can be carried out independently by the training participants, because only what you have already done can be used in the event of danger. The case includes detailed instructions and all relevant utensils for the implementation of 15 different experiments. Chemicals such as baking powder, diesel, petrol, steel wool, et cetera are not included. These must be supplied by the user via the usual procurement paths.

Compliance information

WEEE obligation..... No

Collecting tray mini

Order-No.: 01604



The collecting tub prevents the spread of downcoming liquids. In the event of contamination, the collecting tub can be easily removed and cleaned from the Smartrack®.

Compliance information

WEEE obligation..... No

Collecting tray

Order-No.: 01603



The collecting tub prevents the spread of downcoming liquids. In the event of contamination, the collecting tub can be easily removed and cleaned from the Smartrack®.

Compliance information

WEEE obligation..... No

Storage box insert - for sample vials and test tubes

Order-No.: 01161



Type: 20 x Ø 3 - 10 mm

The module can be inserted precisely into the storage box and enables easy storage and fixing of sample vials and test tubes. It is equipped with a spring mechanism that gives sample bottles and test tubes of different diameters.

Compliance information

WEEE obligation..... No

Storage box insert - for sample vials and test tubes

Order-No.: 01160



Type: 12 x Ø 10 - 20 mm

The module can be inserted precisely into the storage box and enables easy storage and fixing of sample vials and test tubes. It is equipped with a spring mechanism that gives sample bottles and test tubes of different diameters.

Compliance information

WEEE obligation..... No

Storage box insert - for sample vials and test tubes

Order-No.: 01159



Type: 6 x Ø 20 - 30 mm

The module can be inserted precisely into the storage box and enables easy storage and fixing of sample vials and test tubes. It is equipped with a spring mechanism that gives sample bottles and test tubes of different diameters.

Compliance information

WEEE obligation..... No

Storage box insert - for small accessories

Order-No.: 01148



The module can be inserted precisely into the storage box and enables easy storage and fixing of small accessories (e.g. pens, spatula or tweezers).

Compliance information

WEEE obligation..... No

Storage box

Order-No.: 01142



The module is the basis for storing and fixing diverse smaller laboratory equipment. The storage box has a capacity of a maximum of three shelter box inserts that can be arranged and combined individually. Alternatively, the storage box can also be used without shifting box inserts.

Compliance information

WEEE obligation..... No

Storage box insert - for microreaction tubes

Order-No.: 01130



Type: 20 x 2 ml

The module can be inserted precisely into the shelter box and enables easy storage and fixing of microreaction vessels.

Compliance information

WEEE obligation..... No

Storage box insert - for centrifuge tubes

Order-No.: 01129



Type: 6 x 50 ml

The module can be inserted precisely into the storage box and enables easy storage and fixing of centrifuge tube.

Compliance information

WEEE obligation..... No

Storage box insert - for centrifuge tubes

Order-No.: 01128



Type: 12 x 15 ml

The module can be inserted precisely into the storage box and enables easy storage and fixing of centrifuge tube.

Compliance information

WEEE obligation..... No

Storage box insert - for microreaction tubes

Order-No.: 01104



Type: 32 x 0.5 ml

The module can be inserted precisely into the shelter box and enables easy storage and fixing of microreaction vessels.

Compliance information

WEEE obligation..... No

Storage box insert - for centrifuge tubes

Order-No.: 01102



Type: 6 x 15 ml, 4 x 50 ml

The module can be inserted precisely into the storage box and enables easy storage and fixing of centrifuge tube.

Compliance information

WEEE obligation..... No

Additional payment Stock Alternative supplier

Order-No.: 4711SON



Zirconium crucible 35 ml, straight

Order-No.: 11975027



Zirconium tie 35 ml, straight

Centrifuge Super Vario N

Order-No.: 3680



Multi-purpose centrifuge for the dairy industry

Connection value: 230 V/50 ... 60 Hz/1200 VA

Empty weight: 26 kg

Total height: 460 mm with lid

Full height: 370 mm

Speed range: 600 to 1130 UPM **

Tempered

Cellulose swab pure cellulose, 500 pcs/pkg

Order-No.: 47111432123



40x50mm, sterile

Wort broth for microbiology

Order-No.: 1054490500



PU = 500 g

Wort broth for microbiology

Order-No.: 4711WB500



PU = 500 g

Cotton swabs, plastic, sterile

Order-No.: XL561



Material: PP, cotton head made of cotton. Sterile. Eto gas sterilized, individually packed. Medical product class is. Ø head 5 mm. Ø 2.5 mm. Length 150 mm

PU = 100 pieces

Interchangeable membrane heads MC/NH500/TC for ammonium measuring chain

Order-No.: 821212



PU = 5000G

Wort broth for microbiology

Order-No.: 1054495000



PU = 500g

Wort agar for microbiology

Order-No.: 1054480500



Würzebouillon for microbiology

PU = 5000 g

Wort broth for microbiology

Order-No.: 4711WB5000



Strength 102µm, 380 x 508 mm, PU = 100 pieces

Whirl-Pak bag with label, 5441 ml

Order-No.: B01447WA



WH bottle 500ml with red lock PE-HD; Height 136mm; Neckdruckenmeter 58mm;

1 PU = 1 cardboard á 140 doses (1 palette = 12 PU)

WH bottle 500 ml with red cap PE-HD

Order-No.: 200011



Blank, floor and top floor cleared on the outside,

Dimensions: 292 x 257 x 450 mm, Ø 57 mm perforation, weight approx. 1,792 g, inadequate

3A1/Y110, with welded handle

Closure available separately (Article 4711-2004)

Tinplate canister, 30 l

Order-No.: 47115712



Crystal violet neutral red-galle-lactose agar

PU = 500 g

VRB agar for microbiology

Order-No.: 1014060500



For microbiology

PU = 500 g

VRBD agar according to Mossel

Order-No.: 1102750500



Violet Red Bile Dextrose Agar (VRBD -Agar) - EUR. Pharm.

PU = 500g

VRBD agar 500g

Order-No.: 4711VRBD0500



incl. 500 ml bottle, vacusip hand operator, 1- and 8-channel stainless steel suction suction nadel, 1-channel tip adapter f

VACUSIP AKKU aspiration system

Order-No.: 159000



incl. 500 ml bottle, vacusip hand operator, 1- and 8-channel stainless steel suction suction nadel, 1-channel tip adapter for single tips, 1.2 m silicone hose, filter and power supply

VACUSIP suction system

Order-No.: 159010



With sensor, including 4-liter PP bottle, 1 lid with autoclavable PVDF quick couplings and F

VACUSAFE comfort plus

Order-No.: 158310



With glass bottle and sensor, including 3-liter glass bottle, lid with autoclavable PVDF quick couplings and F

VACUSAFE comfort extra

Order-No.: 158300



incl. 4-liter PP bottle, 1 lid with/without plug connections, hose set with protect filter, power supply. Including Vacuboy hand operator, 1-channel stainless steel adapter 40 mm, 1-channel tip adapter with drop, rubber adapter f

VACUSAFE comfort

Order-No.: 158320



Consisting of Vacuboy hand operator, 1-channel and 8-channel stainless steel adapter 40 mm, 1-channel tip adapter f

VACUBOY hand operator set

Order-No.: 155500



including 1-channel stainless steel adapter 40 mm, 1-channel TIP adapter f

VACUBOY hand operator

Order-No.: 155510



For disinfection of room air and surfaces without a stay. Air, surfaces and work surfaces are used by the effect of UV-C light from germs, viruses, yeasts and

Mold spores freed.

Areas of application: production areas I work

UVtec test strips 5-60 mJ/cm²

Order-No.: UVtec560



Urine monovette 10 ml

102 x 15 mm, luer connection, withdrawal tip

PU = 512 pieces

UV-C lamp 900x107x52

Order-No.: 46AR081000



100 g

Urine monovette 10 ml

Order-No.: 10252



For natural gas

Titriplex IV p.a. PU = 100g

Order-No.: 1084240100



PU = 100 pieces

Teclubrenner, 3-flame, with needle valve

Order-No.: 01301E



Diameter above 47 mm

Diameter at the bottom 40 mm

Height 30 mm

Wall thickness 0.035 inches (0.889 mm)

Surface sponge stick SPONGE-STICK with neutralization buffer (NB)

Order-No.: 10006872



Suitable for 1l and 500ml, thread 28 mm

Zirconium melting pot, 35 ml, conical shape

Order-No.: 1035LF



Medium viscos, 200 cst, phenyl -substituted dimethylsiloxan, silicone oil, density (d) 1.04 g/cm³, boiling point (SDP)> 315 ° C, flame point (FLP)> 295 ° C, melting point (F) -55 ° C

PU = 5 kg

Spray attachment for bottles

Order-No.: 4711SHIBSPK1



Special holder Miele laboratory machine for Mojonnier Röhrchen

Silicone oil PN 200, 5 kg

Order-No.: 41302



305x457 mm

PU = 100 St.

Special design

Order-No.: 4711SonANF



With a rooster hose, chrome-plated, suitable connections for fittings, with interior and outer thread M22 and M24

SUL ASSURE PLUS Steri bag

Order-No.: 830251



Connection: M22 x 1, with hose, chrome -plated, length: 150 mm

Not suitable for presseless devices!

Aerator Cascade SLC AC Eco

Order-No.: 47119190422



Sterikon® plus organic indicator for autoclaving control

PU = 100 ampoules

Aerator Cascade SLC

Order-No.: 47111092097



Sterikon® plus organic indicator for autoclaving control

PU = 15 ampoules

Sterikon® plus bioindicator for autoclaving control

Order-No.: 1102740002



Stabilized formazine reference suspensions to determine clarity/opalescence according to European Pharmacopoeia

<0.1; 3; 6; 18; 30 NTU

Sterikon® plus bioindicator for autoclaving control

Order-No.: 1102740001



6 ready-to-use, stabilized formazine primary standard solutions. The set contains sealed kitchen with the following concentrations:

StablCal reference suspension ampoule set

Order-No.: 2897100



From PTFE: Hydrophobe membrane for the filtration of aggressive chemicals, acids and non -polar solvents. Low flow resistance to gases, therefore also used for air and gas filtration (e.g. for ventilation and ventilation tasks and to steri

STABLCAL turbidity standards Calibration set for 2100AN/AN IS

Order-No.: 2659505



With 10 ml neutralization buffer moistened humidified sponge.

PU = 100 sponges per box

No longer available: successor product 700002003

Syringe filter, non-sterile, Ø 15 mm, 0.2µm

Order-No.: KC941



Bag volume: 1000 ml KSB Bagtm

Connection 1: XL hose connection

Connection 2: hose connection (actuator)

Bag: polypropylene (grandpa 15/100 pp)

O-ring: silicone

PU = 100 pieces per cardboard

Sponge-Stick 3M with neutralization buffer

Order-No.: SSL10NB



Product dimensions? 34.5 x 44.2 x 26 cm; 12.61 kilograms

Energy efficiency class? B

Energy efficiency scale (10 classes)

Single-use Keofitt Sampling Bags

Order-No.: 151000



Plastic cubes, external dimensions 72x72x60mm LXBXH, open at the top and bottom, made of 6mm pp, nature, welded

Inner size 60x60x60mm

including loose stamp, approx. 60x60x6/Ø20x120mm

Minimum acceptance quantity: 2 pieces

Siemens iQ300 Microwave

Order-No.: FF023LMB2



Seesand cleaned with acid and glowed for analytical purposes

PU = 25 kg

Set of whipped cream cubes

Order-No.: FT10458



Screw cover vessel 125ml - sterile

Container PP nature, lid LDPE red

1 PU = 380 pieces (1 palette = 12 PU)

Sea sand cleaned with acid and annealed

Order-No.: 4711SA25



Call attachment for a maximum of 36 butyrometer holders or 18 Babcock holders. Radius of the essay: 260 mm

Screw cap jar 125ml - sterile

Order-No.: 4711TP52001



2 connection hoses of 1.5m, connection R3/4 "

Spin attachment A

Order-No.: 3685



Foam roll 160x7 cm

Hose set for ion exchange cartridge

Order-No.: 340037



Role for thermalinter 57 mm x 14 m

PU á 5 pieces

Foam roller 160x7 cm

Order-No.: 4711SN



Ringer tablets for the production of wrestling solution

PU = 100 tablets

Regeneration 2000 with new resin

Order-No.: 4711890027N



Method: Colorimetric with color card 0.04 - 08 - 0.15 ° D 0.05 - 0.1 - 0.19 ° E 0.7 - 1.4 - 2.7 mg/l Caco?
McOlortest TM

400 provisions

Roll for thermoprinter 57 mm x 14 m

Order-No.: 4711TR5714



Regeneration 2800

RINGER tablets for the preparation of RINGER solution

Order-No.: 1155250001



Regeneration 2000

Residual hardness test, Merck

Order-No.: 1111420001



Regeneration 1500

Regeneration 2800

Order-No.: 4711890028



With neutralization buffer. Ready-to-use, pre-dried sponges for easy and easy to extract environment and food samples.

PU = 100 pieces

Regeneration 2000

Order-No.: 4711890027



For the total number determination, 90 mm ready -made plate

PU = 10 pieces

Attention cooling goods! Store at +4 ° C to +12 ° C.

Regeneration 1500

Order-No.: 4711890026



Caseinpepton-Glucose-Hefexttrakt-Agar

According to ISO 4833, ISO 17410 and FDA-BAM Granucult TM

PU = 5 kg

Minimum Shelf Life: 2024/10/21

** Sale - only as long stocks are enough **

Sample bag with sponge stick 38x76 mm

Order-No.: 700002003



PU = 12 pieces

Plate-Count Agar ready-to-use plates

Order-No.: PO5013A



PU = 20 plates

Attention cooling goods! Store at +2 ° C to +8 ° C.

Plate-count agar for microbiology

Order-No.: 1054635000



Blue, 3-layer, 1 roll = 500 demolitions of 38 cm = 190 meters, roll width: 38 cm, point-to-point embossing

PU = 2 roles

Pump tubing, 3-stop, 0.38 mm ID

Order-No.: 47113SPS038



Sample creator 180 ml, with removable handle, sterile - dipper sampler

1 PU = 100 pieces

16 PU per palette

Plate-Count-Agar ready-to-use plates

Order-No.: 1462690020



9mm closure, transparent, shore A, 1.0mm, slit

PU = 100 pieces

Cleaning cloth rolls Multiclean® plus, 3-ply

Order-No.: 1133801



Transparent, 9mm closure, shore A, 1.0mm (agile quality)

PU = 100 pieces

Sample dipper 180 ml, with removable handle sterile - Dipper sampler

Order-No.: 4711SD5302



Pipetboy Pro, white

PP short threaded cap, with hole; silicone white/PTFE blue, 55°

Order-No.: 1200590884



Pipetboy Pro, Pink

PP short threaded cap, with hole; RedRubber / PTFE beige

Order-No.: 1200591003



Pipetboy Pro, Orange

PIPETBOY pro, white

Order-No.: 156400



Pipetboy Pro, green

PIPETBOY pro, pink

Order-No.: 156401



Pipetboy Pro, blue

PIPETBOY pro, orange

Order-No.: 156403



including wall holder, power supply, sterile filter and battery

PIPETBOY pro, green

Order-No.: 156402



including wall holder, power supply, sterile filter and battery

PIPETBOY pro, blue

Order-No.: 156404



including wall holder, power supply, sterile filter and battery

PIPETBOY acu 2 red

Order-No.: 155016



including wall holder, power supply, sterile filter and battery

PIPETBOY acu 2 purple

Order-No.: 155019



including wall holder, power supply, sterile filter and battery

PIPETBOY acu 2 gr

Order-No.: 155015



including wall holder, power supply, sterile filter and battery

PIPETBOY acu 2 color transparent

Order-No.: 155018



For the potentiometric endpoint detection with aging indicator

PIPETBOY acu 2 classic (gray),

Order-No.: 155000



PH inlet stabbing chain SE104N for measurements in semi-solid media, length: 65 /25 mm, shaft material: plastic (Noryl PPO), ceramic diaphragm, pH range: 2 ... 13, temperature range: -5 ... +80 ° C

Sensor connection: DIN 19262

PIPETBOY acu 2 blue

Order-No.: 155017



Including battery, without a measuring chain

pH gel electrode

Order-No.: 11056842



Petri dish PS, 92 x 16 mm, with ventilation cam, sterile

(extremely non-slip / easily stackable)

1 PU = 480 pieces

1 palette = 20 Ve

pH insertion sensor SE 104 N

Order-No.: 616097104



With ventilation cam, easy to stack

PU = 110 pieces in the box (interior packaging: 10)

1 palette = 32 PU = 3520 pieces

Attention: new packaging unit!

pH meter PORTAMESS 911

Order-No.: 616091101



Perasyic acid test 5 - 10 - 20 - 30 - 50 mg/l Merckoquant

PU = 100 pieces

Manufacturer Fa. Merck

Petri dish PS, 92 x 16 mm with cam sterile

Order-No.: 4711BH93B05



Persidic acid test 100 - 150 - 200 - 250 - 300 - 400 - 500 mg/l Merckoquant

PU = 100 pieces

Manufacturer Fa. Merck

Petri dish PS, 140 x 20 mm, sterile

Order-No.: 4711BP14301



Suitable for Coex bottles C40

PERACETIC ACID TEST 5 - 10 - 20 - 30 - 50 MG/L MERCKOQUANT

Order-No.: 1100840001



Reference substance for gas chromatography

5 ml

PERACETIC ACID TEST 100 - 150 - 200 - 250 - 300 - 400 - 500 MG/L MERCCKOQUANT

Order-No.: 1100010001



Microliter syringe Syr 5µl PTFE RN 0.47 (G26S) AS51

Tamper-evident closure C40

Order-No.: 2645



Microliter syringe Syr 1µl TRN 0.64 (G23S) A70ZDV

n-octane 5 ml

Order-No.: 1097160005



12.9 ms/cm at 25 ° C

PU = 250 ml

Multi-injector 7-fold circular 30G / 0.3 x 4 mm

Order-No.: 47111420101



1.41 ms/cm at 25 ° C

PU = 250 ml

Microliter syringe SYR 5µl PTFE RN 0.47 (G26s) AS51

Order-No.: 3030541



147 µs/cm at 25 ° C

PU = 250 ml

Microliter syringe SYR 1µl TRN 0.64(G23s)a70ZDV

Order-No.: 3030609



with piezo ignition

Flame temperature in C °: 1800

Performance in watts: 1650

Gas consumption: 120 g/h

Dimensions with a cartridge: 25 x 11 x 19 cm

Weight: 500g

Matching cartridge: with butangas C 206 (article number: 3218190011)

Conductivity calibration solution 0.1 N

Order-No.: 01NKC1250



Kikkoman for hygiene control

To investigate the contamination of surfaces

ATP/AMP measurement using Lumitester PD20 and PD30.

PU = 100 swabs / box

Conductivity calibration solution 0.01 N

Order-No.: 001NKC1250



Certified reference material for measuring the electrolytic conductivity ($C = 0.1 \text{ mol/l}$) Certipur®

PU = 500 ml

Conductivity calibration solution 0.001 N

Order-No.: 0001NKC1250



Connection $r \frac{3}{4}$ ", measuring range ($\mu S/cm$) 0 - 30 grades LED display/color, pressure -proof up to 6 bar, max. Operating temperature $40^{\circ} C$, voltage supply 12 V DC by means of supplied plug -in power supply, housing polyoxymethylene

Labo-gaz blowtorch soudogazx2000

Order-No.: 3218190071



Analogous 50 for assembly on Tue 1500-15000

Measurement area: 0.1 - 50? S/cm

Connection: R 3/4 "

LuciPac® Pen (for hygiene control)

Order-No.: 1002671



including potential -free contact, assembly mobile or wall mounting

Connection $\frac{3}{4}$ "

Measuring range (? S/cm) 0.055 - 100

Temperature measuring range (° C) 5 - 40

Limit values freely adjustable

Signal at limit value crossing

LF potassium chloride solution (nominal 12.8 mS/cm)

Order-No.: 1012540500



32 x 11.6 mm, wide opening, clear glass, 1.Hydrol. class

PU = 100 pieces

Conductivity meter, analog, with LED

Order-No.: 340044



PU = 2 kg

Conductivity meter R $\frac{3}{4}$ "- analog, with pointer

Order-No.: 340074



Aqualab, WP4C, size: 40mm x 12mm, volume: 15ml

PU = 500 pieces

Conductivity meter C 900, digital

Order-No.: 340057



Air-permeable, bacterial-filtering plugs, ideal for incubation of plant or bacterial cultures as well as for the sterilization of cultural media. Autoclavable up to +200 ° C. Application: Erlenmeyerkroben Enghals, 200-500 ml / Weithals, 50-100 ml. Pa

Short threaded bottle 1.5 ml

Order-No.: 1200588973



Cooling pool fluid 500 ml

Plate count agar (plate count agar)

Order-No.: 4711247910



Fabrtroas to control Uryxxon Relax and Quantofix devices

PU = 1 piece

Plastic sample beaker without lid

Order-No.: 40107



250 ml in PE bottle

Cellulose culture stopper, type 27, 24x36x54 mm

Order-No.: PC241



5G/tablet

PU = 1000 tablets

Cooling bath liquid 500 ml

Order-No.: 874227169



5G/tablet

PU = 250 tablets

Control strip Quantofix Relax

Order-No.: 930988



With grinding gers NS 29/32, Duran® Borosilicate glass 3.3

-Specials-

Control standard C (-0.512°C)

Order-No.: 874237188



Keypoty-free sugar-free fil-IDF for microbiology

PU = 500 g

Kjeldahl tablets (Missouri catalyst)

Order-No.: 1164691000



250 ml in PE bottle

Kjeldahl tablets (Missouri catalyst)

Order-No.: 1164690250



250 ml in PE bottle

Tipping machine attachment 15 ml

Order-No.: 187715



Installation of ion exchanger cartridges

Germ count agar sugar-free FIL-IDF for microbiology

Order-No.: 1108780500



To determine the Titliplex solutions water hardenemit

PU = 500 St.

Calibration standard B (-0.600°C)

Order-No.: 874237187



Hydrogen-peroxide Refill 5510C

Calibration standard A (-0.408°C)

Order-No.: 874237186



Color scale 1200-12,000 ppm

Installation of ion exchanger cartridges

Order-No.: 4711InstDI



Diameter below: 36 mm, diameter above: 44 mm, height: 40 mm, holes: 1 (11 mm Ø)

** Custom -made. Return/exchange excluded **

Indicator buffer tablets

Order-No.: 1084300500



With half hole /1x5 /s, model south milk, for homogenizing pipette (163253022)

PU = 25 pieces

Hydrogen Peroxide Refill 5510C

Order-No.: KAR5510C



Rubber feet for water bath Eurotherm

Hydrogen Peroxide Comparator (flat)

Order-No.: C5510C



Binder -containing, variety MG 227 , diameter: 47 mm

PU = 100 pieces

Rubber stopper NR, gray, Du=36; Do=44; H=40mm

Order-No.: 243410132B



For use in pastless and fatty samples, hard quality,

80 pieces

Rubber stopper 16x21x25mm

Order-No.: DIRI5112500GSSO



Glass fiber fibers, for use in liquid and fatty samples, soft quality, high suction power

PU = 200 pieces

Rubber feet for Eurotherm water bath

Order-No.: AET012



for 2100 N / an / n is / an IS / TL23XX

Glass microfiber filter, MG 227, 47 mm

Order-No.: 301120047



Readability: 0.01 % or 0.001 g

Leader area max.: 101 g

Reproducibility at 2/5 g: 0.1 % / 0.015 %

Drying temperature: 40 ... 230 ° C

Recommended moisture area: 1 ... 100 %

Dimensions (W x D x H): 192 x 385 x 169 mm

Glass fiber fleece, Ø 90 mm

Order-No.: 6906940



CAS no. [7664-39-3]

UN-1790, ADR 8 (6.1) II

PU = 1L

Glass fiber filter Ø 90 mm

Order-No.: 6906941



To store liquid and firm rehearsals. With wire -inserted closure strips. Sterile.

PU = 100 pieces

GELEX test kit, secondary turbidity standard

Order-No.: 2589200



For Lorica 1.2

Moisture analyzer HC103

Order-No.: 30216101



Filter roll for Lamystems Lorica 1.2

Role = 10 m x 0.45 m

Hydrofluoric acid 40 %, p.a., REAG. ISO, REAG. PH. EUR.

Order-No.: 301031L



PU = 80 pieces

Liquid sample bag, 5440 ml, 380x508 mm

Order-No.: 47115440



Glass, Aø 40 x 94mm long, acid and heat -resistant

Filter cell G4 with fiberplast frame

Order-No.: 4711G4N



With corrosion inhibitor

PU = 5L

Filter mat EU 4 PSB/290 S 20/200

Order-No.: 4711FM148



For disinfection control

1. Page: Caso (tryptic soy) agar
2. Page: Caso Agar with enthusiasts

PU = 10 pieces

Disposable aluminum sample dishes

Order-No.: 105013865



For perasy acid sensors

PU = 100ml

Extraction sleeve SE40GL slotted

Order-No.: 130096



6965542

PU = 80 pieces

Ethylene glycol for frost protection

Order-No.: 037705L



Material: PS, polyester fiber plugs. Sterile. Graduation area 1.0 / 2.0 / 2.1 / 2.2 ml. Ø outside 6.3 mm. Top Ø outside 3.7 mm. With cotton stuff

PU = 1000 pieces (40x25)

Envirocheck Contact DC

Order-No.: 1021470001



5 rollers / pack

Electrolyte type PAA

Order-No.: 1023896



- Glass+nickel-plated for viscosimeter Bostwick

Disposable sample dishes, aluminum, round

Order-No.: 645500140



Programmable peristaltics pump, including silicone hose (ID 4 mm, length 2.5 m), focal piece made of stainless steel, stainless steel intake (length 10 and 35 cm), power supply, hose cuff and autoclavable pipette/hose adapter

Demeter pipettes, PS, 2.2 ml sterile

Order-No.: 4711NA441



Transparent, 91x96 mm, screw cover, sterile

PU = 96 pieces

1 palette = 16 Ve

Printer paper for 2100AN and 2100AN IS

Order-No.: 4709000



Article: Debasafe-agriculture

Closure: Debatape

Material: Transparent LDPE 50? M

Printing color: Matt White, Reflex Blue

Print pages: 1 side

Size: 295 x 395 + 30 mm Bon above

Welding seam: pearl seam

Perforation: Leic

Circular spirit level - 25 mm - Clear brass housing

Order-No.: 956369



DAKKKS calibration E2 - 50 g

DOSE-IT P910

Order-No.: 171000



DAKKKS calibration E2 - 200 g

Can with screw cap, PP, 500 ml

Order-No.: 471110191500



DAKKKS calibration E2 - 100 g

DEBAsAFE safety bags 295 x 395 mm

Order-No.: 400516



DAKKKS calibration E2 - 10 g

DAkkS calibration E2 - 50 g

Order-No.: HIG611507



DAKKKS calibration E2 - 1 G

DAkkS calibration E2 - 200 g

Order-No.: HIG611707



PU = 100 pieces

DAkkS calibration E2 - 100 g

Order-No.: HIG611607



Nominal volume 120 ml, material PE/PA, white colored

PU = 360 pieces

DAkkS calibration E2 - 10 g

Order-No.: HIG611307



Nominal volume 550 ml, material PE/PA, white colored

PU = 160 pieces

DAkkS calibration E2 - 1 g

Order-No.: HIG611007



ready-to-use, chromogeneous plate for the detection of bacillus cereus

PU = 240 plates

CEM special crucible 20 ml

Order-No.: 303040



Colorimetric, DPD, with color card and comparener block with long test glasses 0.020 - 0.050 - 0.075 - 0.10
- 0.20 - 0.30 - 0.40 - 0.55 mg/l Clo? McOlortest TM

Coex bottle type C40/OB, 120 ml

Order-No.: 3841



= Myp-agar for microbiology

PU = 500g

Coex bottle type C40/OB, 550 ml

Order-No.: 3660S1



Contact Plates Envirocheck

PU = 20 pieces

Compact Dry X-BC (*Bacillus cereus*)

Order-No.: 1002971



For the residual ceremony specification

Chlorine dioxide test, Merck

Order-No.: 1187540001



Calcium carbonate 99.95 Suprapur VE = 500g

CEREUS selective agar according to Mossel L (BASIS)

Order-No.: 1052670500



Cube with engraving, optionally 1-100

PU = 10 pieces

CASO-AGAR with inhibitors

Order-No.: 1184080020



Container with folding lock / flip-top container PP, 49 ml, aseptic

Fill line at 33ml

1 PU = 600 pieces

1 palette = 30 Ve

Caldur RG 05

Order-No.: 1021000



Filter element with HEPA party filter

Suction power approx. 150m³/h on the suction hood

To recording and filtration of smoke and gases that normally arise when soldering, gluing with solvents and many other work.

Consisting of:

Fan with speed control and filter element

240V, 170W, 50/60Hz

Weight: approx. 8kg fan and filter element

Absaugarm FX2 Ø75 Original L = 1100mm,

Connection below, suction worm with optimized

Flow behavior and integrated barrier flap,

Table attachment, reduction Ø100mm,

Connection hose Ø100 gray, 3m length

Calcium carbonate 99.95 Suprapur PU = 500g

Order-No.: 1020590500



Transparent XL (72 lit, 134 ° C)

PP bag made of particularly strong film with a wide ground seam, for the disposal of contaminated material.

Volume approx. 72 l, foil strength 50? M, bag size 700 mm x 1,100 mm, autoclavable up to 134 ° C

PU = 5

Butter cut strength cubes

Order-No.: 47113170



Returning to SRM from nist nh? cl in H? O 1000 mg/l nh? Certipur®

PU = 500ml

Container with flip-top closure / Flip-Top Container PP, 49 ml

Order-No.: 4711FL5001



Alu foil 45cm x 100; 21My Super P.S. PU (cardboard) = 4 rolls

Suction set EK 3000 EUR 240V EC

Order-No.: 70507844



Pastor of the tube 101 x 16.5 mm, rod material PS, stall length 85 mm, swab material viscose, sterile, 500 pcs/bags

Autoclave bag transparent XL

Order-No.: 54022



transparent XL (72 lit, 134 °C)

PP-Beutel aus besonders starker Folie mit breiter Bodennaht, zur Entsorgung von kontaminiertem Material.

Volumen ca. 72 l, Folienstärke 50 µm, Beutelgröße 700 mm x 1.100 mm, autoklavierbar bis 134°C

VE = 5

Ammonium standard solution, Merck

Order-No.: 1198120500



rückführbar auf SRM von NIST NH_4Cl in H_2O 1000 mg/l NH_4^+ Certipur®

VE = 500ml

Aluminum foil 45cm x 100; 21my Super P.S.

Order-No.: 471118190



Alu-Folie 45cm x 100; 21my Super P.S. VE (Karton) = 4 Rollen

Swab 101x16,5mm, sterile

Order-No.: 80625



Abstrichtupfer im Röhrchen 101 x 16,5 mm, Stabmaterial PS, Stablänge 85 mm, Tupfermaterial Viskose, steril, 500 Stk./Beutel



По вопросам продаж и поддержки обращайтесь:				
Алматы (727)345-47-04 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волгоград (844)278-03-48 Вологда (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89	Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81	Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73 Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47	Ростов-на-Дону (863)308-18-15 Рязань (4912)46-61-64 Самара (846)206-03-16 Санкт-Петербург (812)309-46-40 Саратов (845)249-38-78 Севастополь (8692)22-31-93 Саранск (8342)22-96-24 Симферополь (3652)67-13-56 Смоленск (4812)29-41-54 Сочи (862)225-72-31 Ставрополь (8652)20-65-13 Сургут (3462)77-98-35 Сыктывкар (8212)25-95-17 Тамбов (4752)50-40-97 Тверь (4822)63-31-35	Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удэ (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93
Россия +7(495)268-04-70	Казахстан +7(727)345-47-04	Беларусь +(375)257-127-884	Узбекистан +998(71)205-18-59	Киргизия +996(312)96-26-47

эл.почта: azp@nt-rt.ru || сайт: <https://labor-pilz.nt-rt.ru/>